


Product Datasheet

Rat DLL4 ELISA Kit (orb381197)

Description	Rat DLL4 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S28-E519
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	74961 MW
Uniprot ID	D3ZHH1
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 10-200 pg/ml.
Expiration Date	Please enquire.



Rat DLL4
PicoKine
ELISA Kit
standard
cur...