


Product Datasheet

Rat Angiopoietin-2 ELISA Kit (orb381201)

Description	Rat Angiopoietin-2 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y19-F496
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	56518 MW
Uniprot ID	035462
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is 4-6 ng/ml.
Expiration Date	Please enquire.



Rat
Angiopoietin-
2 PicoKine
ELISA Kit st...