


Product Datasheet

CRISPLD2 antibody (orb38141)

Description	CRISPLD2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2 protein (40-305AA)
Target	CRISPLD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H0B8
Dilution Range	WB:1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(lane 1),...