

Product Datasheet

HSPA1A/HSPA1B antibody (orb38155)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This HSPA1A/HSPA1B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 503-532 amino acids from human HSPA1A/HSPA1B.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P08107
NCBI	NP_005336.3 , NP_005337.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

