



Product Datasheet

CD184/CXCR4 antibody (orb381896)


Description	CD184/CXCR4 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD184/CXCR4. Please contact us for the exact immunogen sequence. The peptide is available as orb13004.
Target	CD184/CXCR4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Dilution Range	IF/ICC: 1:100-800, IHC-P: 1:100-800
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat testis tissue using a...



IHC-P image of rat testis tissue using a...



Immunohistochemical staining of paraffin-embedded tissue using a...