


Product Datasheet

c-Flip protein (orb381903)

Description	c-Flip protein
Conjugation	Unconjugated
Target	c-Flip
Note	For research use only
Application notes	The MW is 55 kDa, and the file stated also 56.5 kDa including his tag.
Protein Sequence	MGSHHHHHHG HMSAEVIHQV EEALDTDEKE MLLFLCRDVA IDVVPPNVRD LLDILRERGK LSVGDLAELL YRVRRFDLLK RILKMDRKAV ETHLLRNPHL VSDYRVLMAE IGEDLDKSDV SSLIFLMKDY MGRGKISKEK SFLDLVVELE KLNLVAPDQL DLLEKCLKNI HRIDLKTKIQ KYKQSVQGAG TSYRNVLQAA IQKSLKDPSN NFRLHNGRSK EQRLKEQLGA QQEPVKKSIQ ESEAFLPQSI PEERYKMMSK PLGICLIIDC IGNETELLRD TFTSLGYEVQ KFLHLSMHGI SQILGQFACM PEHRDYDSFV CVLVSRRGGSQ SVYGVDQTHS GLPLHHIRRM FMGDSCPYLA GKPKMFFIQN YVSEGQLED SSLLEVDGPA MKNVEFKAQK RGLCTVHREA DDFWSLCTAD MSLLEQSHSS PSLYLQCLSQ KLRQERKRPL LDLHIELNGY MYDWNSRVSA KEKYYVWLQH TLRKKLILSY T
MW	56.5 kDa
Expiration Date	6 months from date of receipt.



c-Flip protein