


Product Datasheet

POLR2D antibody (orb381966)

Description	Rabbit polyclonal antibody to POLR2D.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length protein of human POLR2D
Target	POLR2D
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9D7M8 , O15514
Entrez	69241 , 5433
Dilution Range	WB: 1:500-2000, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of POLR2D expressi...



Immunofluorescent analysis of POLR2D sta...