

Product Datasheet

SLC16A4 antibody (orb382041)

Description	Rabbit polyclonal antibody to SLC16A4.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length protein of human MCT4
Target	SLC16A4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8R0M8, O15374
Entrez	229699, 9122
Dilution Range	WB: 1:1000-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.

