

Product Datasheet

SMARCAD1 antibody (orb382045)

Description Rabbit polyclonal antibody to SMARCAD1.

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IF, WB

Immunogen Recombinant full length protein of human SMARCAD1

Target SMARCAD1

Preservatives Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Clonality Polyclonal


Source Rabbit

Uniprot ID [Q9H4L7](#), [Q04692](#), [D3Z9Z9](#)


Entrez [312398](#), [56916](#), [13990](#)

Dilution Range WB: 1:500-2000, IF/ICC: 1:50-100, IP: 1:50-100

Expiration Date 12 months from date of receipt.



Western blot analysis of SMARCAD1 expres...



Immunofluorescent analysis of SMARCAD1 s...