



Product Datasheet

TRPV5 antibody (orb382093)

Description	Rabbit polyclonal antibody to TRPV5.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant full length protein of human TRPV5
Target	TRPV5
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9NQA5 , Q9JIP0 , P69744
Entrez	116469 , 56302 , 194352
Dilution Range	WB: 1:500-2000, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of TRPV5 expression...



Immunofluorescent analysis of TRPV5 staining...