

Product Datasheet

Primate PCSK9 protein (orb3821)

Description	Recombinant of primate PCSK9 protein (C-6His)
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Gln31-Gln152&Ser153-Gln692
Source	Human Cells
Uniprot ID	A8T662
Expiration Date	6 months from date of receipt.