




Product Datasheet

EphA4 antibody (orb38210)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	EPHA4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	109860
Uniprot ID	P54764
NCBI	NP_004429.1
Dilution Range	WB: 1:1000, IF/ICC: 1:100, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of CEM cells us...



Immunofluorescent analysis of human stem c...