

## Product Datasheet

### B7-H6 antibody (orb382134)

## Description

B7-H6 antibody

### Species/Host

Rabbit

### Reactivity

Guinea pig, Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

KLH conjugated synthetic peptide derived from human B7-H6. Please contact us for the exact immunogen sequence. The peptide is available as orb13013.

### Target

B7-H6

### Form/Appearance

10 mM PBS, 0.02% sodium azide

### Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

### Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Polyclonal antibodies are purified by peptide affinity chromatography

### MW

51 kDa

### Uniprot ID

[Q68D85](#)

### NCBI

[001189368](#), [001202439](#)

### Entrez


[374383](#)

### Dilution Range


IHC-P: 1:50-300, WB: 1:200-1000

### Expiration Date


12 months from date of receipt.



IHC-P image of guinea pig brain tissue u...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of human br...