



Product Datasheet

TGN46 antibody (orb382153)


Description	TGN46 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TGN46. Please contact us for the exact immunogen sequence. The peptide is available as orb13033.
Target	TGN46
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	51 kDa
Dilution Range	WB: 1:200-2000, IHC-P: 1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung (lan...



Immunohistochemical staining of paraffin...



IHC-P staining of pig lung tissue using ...