

Product Datasheet

ANO9 peptide (orb382173)

Description	The peptide is used to block anti-ANO9 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ANO9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382458 (ANO9 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P86044, A1A5B4
Entrez	71345, 338440
Expiration Date	12 months from date of receipt.