

Product Datasheet

SRSF5 peptide (orb382177)

Description The peptide is used to block anti-SRSF5 antibody reactivity

Reactivity Human, Mouse, Primate, Rat

Tested Applications Blocking

Target SRSF5

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes This peptide was used as immunogen to produce orb382462 (SRSF5 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID [Q13243](#), [Q09167](#), [O35326](#)

Entrez [20384](#), [6430](#), [29667](#)

Expiration Date 12 months from date of receipt.