

Product Datasheet

ZKSCAN5 peptide (orb382185)

Description	The peptide is used to block anti-ZKSCAN5 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ZKSCAN5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382470 (ZKSCAN5 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z1D8, Q9Y2L8
Entrez	23660, 22757
Expiration Date	12 months from date of receipt.