

Product Datasheet

GTPBP2 peptide (orb382193)

Description	The peptide is used to block anti-GTPBP2 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	GTPBP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382478 (GTPBP2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q3UJK4, Q9BX10
Entrez	54676, 56055
Expiration Date	12 months from date of receipt.