

Product Datasheet

SLC30A2 peptide (orb382210)

Description	The peptide is used to block anti-SLC30A2 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SLC30A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382495 (SLC30A2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q2HJ10 , Q62941 , Q9BRI3
Entrez	230810 , 25362 , 7780
Expiration Date	12 months from date of receipt.