

Product Datasheet

XRCC6 peptide (orb382224)

Description	The peptide is used to block anti-XRCC6 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	XRCC6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382509 (XRCC6 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23475, P12956
Entrez	14375, 2547
Expiration Date	12 months from date of receipt.