

Product Datasheet

MYH1 peptide (orb382227)

Description	The peptide is used to block anti-MYH1 antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	MYH1; MYH2; MYH3; PTPRJ; MYH6; MYH7; MYH8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382512 (MYH1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UKX2 , Q9Y623 , Q5SX39 , P12882 , P12847 , P11055 , P13533 , P13541 , P13535 , Q29RW1 , P12883 , Q5SX40
Entrez	4619 , 4621 , 17884 , 360543 , 17879 , 4624 , 17883 , 5795 , 24583 , 4620 , 4625 , 4626
Expiration Date	12 months from date of receipt.