

Product Datasheet

APOA5 peptide (orb382245)

Description	The peptide is used to block anti-APOA5 antibody reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	APOA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382528 (APOA5 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9QUH3, Q6Q788, Q8C7G5
Entrez	140638, 116519, 66113
Expiration Date	12 months from date of receipt.