

Product Datasheet

ARL2BP peptide (orb382258)

Description	The peptide is used to block anti-ARL2BP antibody reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ARL2BP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382541 (ARL2BP antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D385 , Q4V8C5 , Q9Y2Y0
Entrez	23568 , 107566 , 498910
Expiration Date	12 months from date of receipt.