

Product Datasheet

TBX2 peptide (orb382269)

Description	The peptide is used to block anti-TBX2 antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	TBX2; TBX3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382552 (TBX2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13207 , Q7TST9 , Q60707 , O15119 , P70324
Entrez	6926 , 353305 , 21386 , 21385
Expiration Date	12 months from date of receipt.