

Product Datasheet

GDI2 peptide (orb382290)

Description	The peptide is used to block anti-GDI2 antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	GDI2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382570 (GDI2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61598, P50395, P50399
Entrez	14569, 29662, 2665
Expiration Date	12 months from date of receipt.