

Product Datasheet

IDH3A peptide (orb382291)

Description	The peptide is used to block anti-IDH3A antibody reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	IDH3A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382571 (IDH3A antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9D6R2 , P50213 , Q99NA5
Entrez	67834 , 3419 , 114096
Expiration Date	12 months from date of receipt.