

Product Datasheet

TYRO3 peptide (orb382293)

Description	The peptide is used to block anti-TYRO3 antibody reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	TYRO3; MERTK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382573 (TYRO3 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P55146 , P55144 , Q06418 , P57097 , Q12866 , Q60805
Entrez	22174 , 25232 , 17289 , 7301 , 65037 , 10461
Expiration Date	12 months from date of receipt.