

Product Datasheet

MERTK peptide (orb382294)

Description	The peptide is used to block anti-MERTK antibody reactivity
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MERTK; TYRO3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382574 (MERTK antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P55146, P55144, Q06418, P57097, Q12866, Q60805
Entrez	22174, 25232, 17289, 7301, 65037, 10461
Expiration Date	12 months from date of receipt.