

Product Datasheet

WIPI1 peptide (orb382303)

Description	The peptide is used to block anti-WIPI1 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	WIPI1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382583 (WIPI1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8R3E3, Q5MNZ9
Entrez	52639, 55062
Expiration Date	12 months from date of receipt.