

Product Datasheet

HIST1H2BC peptide (orb382314)

Description	The peptide is used to block anti-HIST1H2BC antibody reactivity
Reactivity	Human, Mouse, Primate
Tested Applications	Blocking
Target	HIST1H2BC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382593 (HIST1H2BC antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62807, P70696, P33778, Q6ZWY9, Q64475, Q96A08
Entrez	8346, 319177, 8339, 8343, 8347, 319179, 3017, 319178, 3018, 8344, 255626
Expiration Date	12 months from date of receipt.