

Product Datasheet

LGALS3BP peptide (orb382331)

Description	The peptide is used to block anti-LGALS3BP antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	LGALS3BP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382610 (LGALS3BP antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q07797 , Q08380
Entrez	19039 , 3959
Expiration Date	12 months from date of receipt.