

Product Datasheet

CSNK1A1 peptide (orb382369)

Description	The peptide is used to block anti-CSNK1A1 antibody reactivity
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	CSNK1A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382647 (CSNK1A1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BK63 , P97633 , P48729
Entrez	93687 , 113927 , 1452
Expiration Date	12 months from date of receipt.