

Product Datasheet

KRT14 peptide (orb382375)

Description	The peptide is used to block anti-KRT14 antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	KRT14; KRT17
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382653 (KRT14 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P02533, Q6IFV1, Q6IFU8, Q61781, Q9QWL7, Q04695
Entrez	16664, 16667, 287702, 3872, 3861, 287701
Expiration Date	12 months from date of receipt.