

Product Datasheet

ACTA2 peptide (orb382385)

Description	ACTA2 peptide
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	ACTA2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382726 (ACTA2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62736, P62738, P68032, P63261, P63267, P62737, P68133
Entrez	11475, 81633, 59
Expiration Date	12 months from date of receipt.