

Product Datasheet

TGFB1I1 peptide (orb382390)

Description	The peptide is used to block anti-TGFB1I1 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	TGFB1I1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382667 (TGFB1I1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O43294, Q62219
Entrez	21804, 7041
Expiration Date	12 months from date of receipt.