

Product Datasheet

C8B peptide (orb382397)

Description	The peptide is used to block anti-C8B antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	C8B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382674 (C8B antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P07358, Q8BH35
Entrez	732, 110382
Expiration Date	12 months from date of receipt.