

Product Datasheet

MYOZ2 peptide (orb382402)

Description	The peptide is used to block anti-MYOZ2 antibody reactivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	MYOZ2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382679 (MYOZ2 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NPC6, Q9JJW5
Entrez	59006, 51778
Expiration Date	12 months from date of receipt.