

Product Datasheet

Neurokinin B receptor Peptides (orb382415)

Description	Neurokinin B receptor Peptides
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TACR3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382692 (TACR3 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P16177, P47937, P29371
Entrez	6870, 24808, 21338
Expiration Date	12 months from date of receipt.