

Product Datasheet

Collagen I peptide (orb382424)

Description	Collagen I peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	COL1A1; COL1A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382699 (COL1A1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P02466, P08123, P02454, P02452, P11087, Q01149
Entrez	12843, 1277, 1278, 29393, 84352, 12842
Expiration Date	12 months from date of receipt.