

Product Datasheet

MAP1LC3B peptide (orb382441)

Description	The peptide is used to block anti-MAP1LC3B antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAP1LC3B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382715 (MAP1LC3B antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CQV6 , Q9GZQ8 , Q62625
Entrez	67443 , 81631 , 64862
Expiration Date	12 months from date of receipt.