


Product Datasheet

LRTM1 antibody (orb38248)


Description	LRTM1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Leucine-rich repeat and transmembrane domain-containing protein 1 protein (1-212AA)
Target	LRTM1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HBL6
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:10-100
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(lane 1), H...



Immunohistochemical staining of human pa...



Immunofluorescence analysis of Hela cell...