


Product Datasheet

MYH1 antibody (orb382512)

Description	Rabbit polyclonal antibody to MYH1
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K505 of human Myosin-pan protein. The exact sequence is proprietary.
Target	MYH1; MYH2; MYH3; PTPRJ; MYH6; MYH7; MYH8
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P12847 , Q5SX39 , Q9UKX2 , P12883 , P12882 , P13533 , Q9Y623 , P11055 , P13535 , Q5SX40 , Q29RW1 , P13541
Entrez	17879 , 4625 , 24583 , 360543 , 17883 , 4619 , 17884 , 4624 , 4626 , 5795 , 4620 , 4621
Dilution Range	WB: 1:500:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Myosin-pan (Ace...)