


Product Datasheet

TUBB2A antibody (orb382514)

Description	Rabbit polyclonal antibody to TUBB2A
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic acetylated peptide corresponding to residues surrounding K379 of human Beta-tubulin protein. The exact sequence is proprietary.
Target	TUBB2A; TUBB2B; TUBB3; TUBB4A; TUBB4B; TUBB; TUBB6
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9D6F9 , P68371 , Q4QRB4 , P07437 , Q7TMM9 , P85108 , Q13509 , Q13885 , Q9ERD7 , Q9BVA1 , Q3KRE8 , Q9CWF2 , P04350 , Q9BUF5
Entrez	10383 , 7280 , 84617 , 203068 , 22153 , 22151 , 10381 , 22152 , 347733 , 10382 , 246118 , 291081 , 73710 , 498736
Dilution Range	WB: 1:500:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Beta-tubulin (A...