



Product Datasheet

ASF1A antibody (orb382526)

Description	Rabbit polyclonal antibody to ASF1A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ASF1A. The exact sequence is proprietary.
Target	ASF1A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9Y294 , Q9CQE6
Entrez	25842 , 66403
Dilution Range	WB: 1:500:2000, IF/ICC: 1:50:100
Expiration Date	12 months from date of receipt.



Western blot analysis of ASF1A expression...



Immunofluorescent analysis of ASF1A staining...