



Product Datasheet

TYRO3 antibody (orb382573)

Description	Rabbit polyclonal antibody to TYRO3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y749/681 of human MERTK/TYRO3 protein. The exact sequence is proprietary.
Target	MERTK; TYRO3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P57097 , P55146 , P55144 , Q12866 , Q60805 , Q06418
Entrez	65037 , 7301 , 17289 , 22174 , 25232 , 10461
Dilution Range	WB: 1:500:2000, IHC-P: 1:50:200
Expiration Date	12 months from date of receipt.



Western blot analysis of MERTK/TYRO3 (Phosphorylated)



Immunohistochemical analysis of MERTK/TYRO3