

Product Datasheet

KRT2 antibody (orb382594)

Description	Rabbit polyclonal antibody to KRT2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Cytokeratin-pan. The exact sequence is proprietary.
Target	KRT2; KRT76; KRT3; KRT5; KRT6A; KRT6B; KRT6C; KRT71; KRT72; KRT73; KRT74; KRT75; KRT79; KRT8; KRT84
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P04259, P08729, Q86Y46, Q3SY84, Q14CN4, P48668, Q7RTS7, O95678, Q5XKE5, P02538, P12035, Q922U2, Q9NSB2, Q6P6Q2, Q01546, P35908, Q3TTY5, Q3UV17, Q6IG02, P13647, P05787
Entrez	338785, 112802, 3856, 9119, 319101, 3890, 3849, 406228, 16681, 3850, 140807, 3853, 3852, 110308, 77055, 369017, 51350, 286887, 3854, 121391
Dilution Range	WB: 1:500:2000
Expiration Date	12 months from date of receipt.

