



## Product Datasheet

### UBA52 antibody (orb382596)

|                            |                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to UBA52                                                                                               |
| <b>Species/Host</b>        | Rabbit                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rat, Sheep                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                      |
| <b>Tested Applications</b> | IH, WB                                                                                                                            |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Ub. The exact sequence is proprietary. |
| <b>Target</b>              | UBA52; RPS27A; UBB; UBC                                                                                                           |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                         |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                         |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                  |
| <b>Note</b>                | For research use only                                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                                                        |
| <b>Source</b>              | Rabbit                                                                                                                            |
| <b>Uniprot ID</b>          | <b>P62987, P0CG47, P62979, P0CG48</b>                                                                                             |
| <b>Entrez</b>              | <b>22186, 6233, 192255, 7311, 7314, 7316, 50522, 22190, 64156, 100912032, 22187, 78294</b>                                        |
| <b>Dilution Range</b>      | WB: 1:500:2000, IHC-P: 1:50:200                                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                   |



Western blot analysis of Ub expression i...



Immunohistochemical analysis of Ub stain...