




Product Datasheet

HER2 antibody (orb38262)

Description	Rabbit polyclonal antibody to HER2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ERBB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	137910
Uniprot ID	P04626
NCBI	NP_004439.2 , NP_001005862.1
Dilution Range	WB: 1:1000, FACS: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of ErbB2 antibody ...



Flow cytometric analysis of HER2/ErbB2 a...