


Product Datasheet

TLE1 antibody (orb382641)

Description	Rabbit polyclonal antibody to TLE1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human TLE1/2/3/4. The exact sequence is proprietary.
Target	TLE1; TLE2; TLE3; TLE4
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q62441 , Q04725 , Q62440 , Q9WVB2 , Q04727 , Q08122 , Q07141 , Q04726 , Q04724 , Q9JIT3
Entrez	21888 , 7088 , 7091 , 21887 , 25565 , 21885 , 21886 , 84424 , 7090 , 7089
Dilution Range	WB: 1:500:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of TLE1/2/3/4 expr...