



Product Datasheet

YARS antibody (orb382663)

Description	Rabbit polyclonal antibody to YARS
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Tyrosyl-tRNA Synthetase. The exact sequence is proprietary.
Target	YARS
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q4KM49 , P54577 , Q91WQ3
Entrez	107271 , 313047 , 8565
Dilution Range	WB: 1:500:2000, IHC-P: 1:50:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Tyrosyl-tRNA Synthetase (YARS) in two lanes, A and B.



Immunohistochemical analysis of Tyrosyl-tRNA Synthetase (YARS) in tissue sections.